

**LC.** 3D HPLC system consists of a pumping system, injection valve, column, trifunctional detector, and necessary tubing and connectors. Simultaneous detection in UV, conductivity, and fluorescence modes is possible. HPLC Technology 409

**STM.** Tunnelscope 2400 automatic scanning tunneling microscope scans areas up to  $10 \times 10 \mu\text{m}$  with a vertical resolution better than 0.01 nm and typical lateral resolution of 0.01 nm. Coarse and fine positioning of the tip are electronically controlled. Struers 411

**Cryogenics.** Model 9000 Controller is a cryogenic freezing system that allows full control of all freezing parameters through screen-prompted commands. The system can generate graphs showing the temperature of the sample and the freezing chamber as the specimen freezes. Gordinier Electronics 426

**Spectrophotometers.** DU series 7000 diode array spectrophotometers feature microanalysis capabilities; applications software; and full spectrum quantitation, in which sample information from the entire scanned spectrum

is used to quantitate the concentration of individual components in a mixture. Beckman Instruments 427

## Software

**Surface analysis.** SIMS software, which operates under a Microsoft Windows 3.0 environment, features Spearman's rank correlation; addition, subtraction, and overlay of SIMS data; and a plotting routine that allows multiple plots on a page. Sub Monolayer Science 431

**Data acquisition.** Easyest is an icon-driven data acquisition and analysis software package designed for IBM PCs and compatible computers. Features include simplified parameter setting and multiple-channel monitoring. Keithley Asyst 432

## Manufacturers' Literature

**SEC.** Brochure describes high-efficiency columns for size-exclusion chromatography. Resin-, silica- and styrene-divinyl benzene-based columns are discussed. Supelco 429

**Environmental analysis.** Three application notes discuss LC, GC, and GC/MS methods for the analysis of drinking water and wastewater samples. Perkin-Elmer 430

**Dioxin.** Brochure highlights the DB-Dioxin GC column, which separates 2,3,7,8-TCDD and 2,3,7,8-TCDF from all other tetra-substituted dioxins and furans. J&W Scientific 434


**LC.** *Mac-Mod Forum* contains information on changes in column selectivity and temperature effects in HPLC. A column for the separation of enantiometric compounds is also described. Mac-Mod Analytical 435

**Electrochemistry.** Application note AC-5 provides a survey of applications of electrochemical impedance in the evaluation of coatings on metals. EG&G Princeton Applied Research 436

**Spectrophotometer.** Brochure highlights the Model F-3010 fluorescence spectrophotometer, which features a horizontal optical beam configuration and monochromatic ratioing system. Hitachi 437


# Malvern - the complete answer to particle characterization

Malvern offer the most comprehensive and advanced range of particle sizing and zeta potential equipment available anywhere.



Zeta potential and submicron particle sizing - the Malvern ZetaSizer 3 offers the ultimate in PCS particle sizing using unique multi-angle analysis for reliable, detailed size distributions from 3 nanometres to 3 microns. Its electrophoresis capabilities extend from strong electrolyte media to non-polar liquids such as hydrocarbons with full zeta potential distributions measured in seconds. For characterizing colloidal dispersions, it has no equal.

Direct dry powder analysis with MasterSizer - already the world's most popular laser particle sizer, its new Dry Powder Feeder makes it even quicker and more versatile. Overall size range is 0.1-600 microns and MasterSizer is the only one to offer genuine Mie analysis.



The Malvern API AeroSizer uses aerodynamic techniques combined with advanced laser technology to give extremely high resolution particle size analysis of powders from 0.5 to 200 microns.



Malvern Instruments Inc.  
10 Southville Rd, Southborough, MA 01772, USA  
Tel: (508) 480 0200 Fax: (508) 460 9692

Malvern Instruments Limited  
Spring Lane South, Malvern, Worcs, WR14 1AQ, UK  
Telephone: (0684) 892456 Telex: 339679 Fax: (0684) 892789

CIRCLE 92 ON READER SERVICE CARD