

AUGUST 10-14

Molecular Theory and Simulation of Adsorption

P. A. Monson, *discussion leader*

Computer Simulation of Gases Sorbed in Realistic Model Pores. W. A. Steele

Application of Density Functional Theories to Adsorption of Simple Fluid Mixtures. M. L. Rosinberg

Molecular Modeling of Adsorption and Transport in Zeolite Molecular Sieves. D. N. Theodorou

Thermodynamic Theories of Adsorption of Mixtures Based on Molecular Simulations. A. L. Myers

Adsorption Separation

R. L. Albright, *discussion leader*

New Adsorbents for Air Purification. C. Chang, N. N. Li

Pressure Swing Adsorption: The Next Distillation? K. S. Knaebel

Kinetically Controlled PSA Separations: Comparison with Membrane Separations. D. M. Ruthven

Membrane Separation

B. Bikson, *discussion leader*

Olefin Separations in Ion-Exchange Membranes. R. D. Noble

Nanoparticulates, Nanocavities, and Nanospace Separation in Ceramic Membranes. M. A. Anderson

Process Models for Control of Separations Processes

J. R. Fair, *discussion leader*

Nonlinear Wave Theory for Dynamics of High-Purity Distillation Columns. Y.-L. Hwang

Dynamic Behavior and Multiple Steady States in Azeotropic Distillation Sequences. M. F. Doherty

Using Tray-to-Tray Models for the Control of Distillation Columns. J. B. Riggs

Use of Steady-State and Dynamic Process Models in the Development of Process Control Strategies: An Industrial Perspective. J. J. Downs

Bioseparation

J. C. Baker, *discussion leader*

Critical Displacer Properties in Protein Displacement Chromatography. S. M. Cramer

Competitive Protein Adsorption in Large-Scale Affinity Chromatography. N.-H. L. Wang

Separation and Focusing Biochemicals in Countercurrent Magnetically Stabilized, Fluidized Bed Adsorption Columns. M. A. Burns

Environmental Applications

New Separation Opportunities in

Environmental Applications. G. E. Keller

New Commercial Processes

G. L. Hubred, *discussion leader*

Solvent Pulping: Promise and Performance. R. Katzen

Coupling the Separation Technologies into BHP's Resource Processes: The Engineering-Economic Dimensions. G. Rigby

Immunoaffinity Process for Production of Activated Protein C, A New Plasma-Derived Antithrombotic Agent. C. L. Orthner

Microwave Digestion Bombs

In two sizes, 23 ml and 45 ml.



Rapid microwave heating can now be used to speed the treatment of inorganic, organic and biological samples in strong acid or alkali media using Parr Microwave Digestion Bombs. These bombs have a microwave-transparent outer body which encloses an inner Teflon sample holder, allowing microwave energy to be used to develop high internal temperatures and pressures without overheating the outer wall. Samples can be dissolved with heating times of one minute or less with complete safety when the bombs are used as directed. Call or write for Bulletin 4700 for complete technical details.



PARR INSTRUMENT COMPANY
211 Fifty-third Street • Moline, IL 61265 • Telex: 270226
Phone: 309-762-7716 • 800-872-7720 • Fax: 309-762-9453

CIRCLE 108 ON READER SERVICE CARD

ACS Software

High-quality software programs for the personal computer that meet the standards you expect from the ACS

For information about ACS SOFTWARE, call TOLL FREE (800) 227-5558 or write to:
American Chemical Society,
Marketing Communications,
1155 Sixteenth Street, N.W.
Washington, D.C. 20036

Today's Best Source of Quantitative Numerical Data of Physics and Chemistry JOURNAL OF Physical and Chemical Reference Data

Editor, David R. Lide, Jr.
National Institute of Standards
and Technology

Published BIMONTHLY by the American Chemical Society, the American Institute of Physics, and the National Institute of Standards and Technology, JPCRD provides you with compilations and reviews produced under the National Standard Reference Data System of the National Institute of Standards and Technology.

JPCRD contains recommended values, uncertainty limits, critical commentary on methods of measurement, and full references to the original papers.

Call Toll Free 1-800-333-9511 and charge your order! Outside the U.S. call (614) 447-3776, or write:

American Chemical Society,
Member and Subscriber Services
P.O. Box 3337
Columbus, OH 43210 U.S.A.